



QUICK REFERENCE GUIDE

PACKAGING, CONFIGURATION & STORAGE

KIT CONTAINS:

Before you start, review the contents of the kit first and read the instructions carefully.

- Test Cassette – 20 each – An absorbent membrane cassette individually wrapped in foil pouch, containing a plastic pipette for blood sample.
- Dropper Bottle – 6 mL – containing Buffer diluent.
- Spring loaded lancets – 20 each.
- Positive Control – 1 each – Positive TSH serum for quality control testing of the system (See Quality Control Section).
- Negative Control – 1 each – Negative TSH serum for quality control testing of the system (see Quality Control Section).

ADDITIONAL COMPONENTS REQUIRED (not provided in kit)

- Timer
- Alcohol wipe
- External Controls (optional)

KIT STORAGE

- The test kit may be stored at room temperature (15 – 30°C; 60 – 86°F); **do not freeze.**
- Do not use the test cassette after the date printed on the foil pouch.
- Keep away from moisture, heat, or direct sunlight